




## YC-400L

# Best Seller

## 2~8°C Pharmacy Refrigerator

Seven Sensors  
 Temperature Uniformity  
 Ensuring Safety Storage



Temp. Uniformity



Pulling Down to 5°C



Temp. Recovery After Door Opened for 1 Min

*\*Typical freezer data based on internal testing with freezer set point at 5°C and ambient temperature at 22°C.*

## Features



**Temperature Control System**  
 Display temp. in 0.1°C increments



**Double-pane Glass Door**  
 Automatic electric heating function



**Air Cooling System**



**CE Certificate**




**Foaming Door**



**Alarms System**  
 Multiple alarm functions

## Structure





| 2~8°C Pharmacy Refrigerator   |                                                                                  |                  |
|-------------------------------|----------------------------------------------------------------------------------|------------------|
| Model                         | YC-400L                                                                          |                  |
| Capacity(L)                   | 400                                                                              |                  |
| Internal Size (W*D*H) mm/Inch | 604×520×1400/23.8*20.5*55.1                                                      |                  |
| External Size (W*D*H) mm/Inch | 700*645*2010/27.6*25.4*79.1                                                      |                  |
| Package Size (W*D*H) mm/Inch  | 765*710*2074/30.1*28.0*81.7                                                      |                  |
| NW/GW (kgs)                   | 116/130                                                                          |                  |
| Performance                   |                                                                                  |                  |
| Temperature Range             | 2~8°C                                                                            |                  |
| Ambient Temperature           | 16~32°C                                                                          |                  |
| Cooling Performance           | 5°C                                                                              |                  |
| Climate Class                 | N                                                                                |                  |
| Controller                    | Microprocessor                                                                   |                  |
| Display                       | Digital Display                                                                  |                  |
| Refrigeration                 |                                                                                  |                  |
| Compressor                    | 1pc                                                                              |                  |
| Cooling Method                | Air Cooling                                                                      |                  |
| Defrost Mode                  | Automatical                                                                      |                  |
| Refrigerant                   | R290                                                                             |                  |
| Insulation Thickness(mm)      | D:55, L/R:48, B:50                                                               |                  |
| Construction                  |                                                                                  |                  |
| External Material             | PCM                                                                              |                  |
| Inner Material                | Stainless Steel                                                                  |                  |
| Shelves                       | 6                                                                                |                  |
| Access Test Port              | 1pc. Ø 23mm                                                                      |                  |
| Casters                       | 4(2 brakes)                                                                      |                  |
| Data Logging/Export           | Record every 10 mins/Max.12 months via USB port                                  |                  |
| Interface                     | Remote Alarm Contact, RS485                                                      |                  |
| Backup Battery                | Yes                                                                              |                  |
| Alarms                        |                                                                                  |                  |
| Temperature                   | High/Low Temperature, High Ambient Temperature, Condenser High Temperature Alarm |                  |
| Electrical                    | Power Failure, Low Battery, Sensor Error                                         |                  |
| System                        | Data Communication Failure, Door Ajar,Data Logger Failure                        |                  |
| Electrical                    |                                                                                  |                  |
| Power Supply(V/Hz)            | 110V/60Hz                                                                        | 220-240V;50/60Hz |
| Rated Current(A)              | 4.74A                                                                            | 2.23A            |

\*Zhongke Meiling reserves the right to change products and specifications without prior notice.